

Table of contents

\\NEUROSPIN

<table border="1"><tr><td>RECHERCHE</td></tr></table>	RECHERCHE
RECHERCHE	

<table border="1"><tr><td>FM</td></tr></table>	FM
FM	

<table border="1"><tr><td>C2P_Protocols_QA</td></tr></table>	C2P_Protocols_QA
C2P_Protocols_QA	

<table border="1"><tr><td>1-epi_QA_Noise_SNR</td></tr></table>	1-epi_QA_Noise_SNR
1-epi_QA_Noise_SNR	

<table border="1"><tr><td>2-epi_QA_Stability_warmup_no_reco</td></tr></table>	2-epi_QA_Stability_warmup_no_reco
2-epi_QA_Stability_warmup_no_reco	

<table border="1"><tr><td>3-epi_QA_Stability</td></tr></table>	3-epi_QA_Stability
3-epi_QA_Stability	

\NEUROSPIN\RECHERCHE\FM\C2P_Protocols_QA\1-epi_QA_Noise_SNR

TA: 10 sec Coil Selection: Auto Voxel Size: 4.0×4.0×4.0 mm³ Acc:: None Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Resolution - Common

Base Resolution	64
Phase Resolution	50 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	None
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	On
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Normalize	Off

Routine

Slice Group	1
Slices	1
Distance Factor	50 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	256 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm
TR	5000.0 ms
TE	8.80 ms
Averages	1
Concatenations	1
AutoAlign	---

Geometry - Common

Slice Group	1
Slices	1
Distance Factor	50 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	256 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm
TR	5000.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Contrast - Common

TR	5000.0 ms
TE	8.80 ms
MTC	Off
Flip Angle	80 deg
Fat-Water Contrast	Standard
Reconstruction	Magnitude

Geometry - AutoAlign

Slice Group	1
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	2
Multiple Series	Off
Delay in TR	0.00 ms

Geometry - Saturation

Special Saturation	None
--------------------	------

Resolution - Common

FOV Read	256 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm

Geometry - Tim Planning Suite

Table Position	H
----------------	---

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.159149 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	5000.0 ms
Concatenations	1

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Excitation	Standard
RF Pulse Type	Fast
Gradient Mode	Fast

Sequence - Part 1

Bandwidth	4882 Hz/Px
Echo Spacing	0.35 ms
Free Echo Spacing	Off
EPI Factor	32

Sequence - Part 2

Introduction	Off
--------------	-----

Sequence - Special

QA Postprocessing	Noise & SNR
Mosaic	On

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\NEUROSPIN\RECHERCHE\FM\C2P_Protocols_QA\2-epi_QA_Stability_warmup_no_reco

TA: 6:44 min Coil Selection: Auto Voxel Size: 3.4×3.4×4.0 mm³ Acc:: None Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Resolution - Common

Base Resolution	64
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	None
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	On
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Normalize	Off

Routine

Slice Group	1
Slices	27
Distance Factor	25 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm
TR	2000.0 ms
TE	26.00 ms
Averages	1
Concatenations	1
AutoAlign	---

Geometry - Common

Slice Group	1
Slices	27
Distance Factor	25 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm
TR	2000.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Contrast - Common

TR	2000.0 ms
TE	26.00 ms
MTC	Off
Flip Angle	90 deg
Fat-Water Contrast	Fat Saturation
Reconstruction	Magnitude

Geometry - AutoAlign

Slice Group	1
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	200
Multiple Series	Off
Delay in TR	0.00 ms

Geometry - Saturation

Special Saturation	None
--------------------	------

Resolution - Common

FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm

Geometry - Tim Planning Suite

Table Position	H
----------------	---

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	220 mm
R >> L	220 mm
F >> H	134 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.159149 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	2000.0 ms
Concatenations	1

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Excitation	Standard
RF Pulse Type	Fast
Gradient Mode	Fast

Sequence - Part 1

Bandwidth	1698 Hz/Px
Echo Spacing	0.65 ms
Free Echo Spacing	Off
EPI Factor	64

Sequence - Part 2

Introduction	Off
--------------	-----

Sequence - Special

QA Postprocessing	None
Mosaic	On

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\NEUROSPIN\RECHERCHE\FM\C2P_Protocols_QA\3-epi_QA_Stability

TA: 6:44 min Coil Selection: Auto Voxel Size: 3.4×3.4×4.0 mm³ Acc:: None Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Resolution - Common

Base Resolution	64
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	None
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	On
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Normalize	Off

Routine

Slice Group	1
Slices	27
Distance Factor	25 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm
TR	2000.0 ms
TE	26.00 ms
Averages	1
Concatenations	1
AutoAlign	---

Geometry - Common

Slice Group	1
Slices	27
Distance Factor	25 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm
TR	2000.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Contrast - Common

TR	2000.0 ms
TE	26.00 ms
MTC	Off
Flip Angle	90 deg
Fat-Water Contrast	Fat Saturation
Reconstruction	Magnitude

Geometry - AutoAlign

Slice Group	1
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	200
Multiple Series	Off
Delay in TR	0.00 ms

Geometry - Saturation

Special Saturation	None
--------------------	------

Resolution - Common

FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm

Geometry - Tim Planning Suite

Table Position	H
----------------	---

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	220 mm
R >> L	220 mm
F >> H	134 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.159149 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	2000.0 ms
Concatenations	1

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Excitation	Standard
RF Pulse Type	Fast
Gradient Mode	Fast

Sequence - Part 1

Bandwidth	1698 Hz/Px
Echo Spacing	0.65 ms
Free Echo Spacing	Off
EPI Factor	64

Sequence - Part 2

Introduction	Off
--------------	-----

Sequence - Special

QA Postprocessing	EPI Stability
Mosaic	On

Sequence - Assistant

SAR Assistant	Off
---------------	-----