

# **QA for pTx coils**

Franck Mauconduit

DRF/Joliot/NeuroSpin & Université Paris-Saclay, CEA, Gif-sur-Yvette, France

*March 3rd 2023*


# Methods

- Method : acquisition of B1+ maps and comparison to a reference B1+ map
- Installation :
  - Predefined phantom : Agar (no fluid movement)  
short T1 : ~650ms
  - Reproducible positioning : holder available + automatic position detection  
with a profile measurement in 3 main directions
- Sequence : based on an XFL acquisition  
(Fast mode = 9 gre + 1 saturated gre)
- Protocol : must be fixed in advanced for comparison and rapid to execute  
5 slices, TR = 2s, XFL fast mode => 20s total acquisition time (for 8 Tx channels)



# Methods

List of protocols used for the QA of a pTx coil



# Reconstruction process


Reference file :

location : %CustomerSeq%/SafetyCoilCheck/

filename : [COILNAME]\_fast\_reference.dat

ex : 8Tx32Rx\_Head\_7T\_fast\_reference.dat

→ Each coil will have its reference file



Measurement file :

Standard reconstruction from our B1 mapping XFL sequence

location : %CustomerSeq%/RECO\_B0B1/

This measurement is compared to the reference file in ICE

Results :

Dicom series :

- [PROTOCOL\_NAME]\_CHECK-PTX-COIL  
phase images -> Visual check of the magnitude and
- [PROTOCOL\_NAME]\_FOLLOWUP-PTX-COIL  
measurement -> Follow-up graphs and validity of the

Followup data saved in %CustomerSeq%/QA/[COILNAME]/


# Visual check

Dicom series [PROTOCOL\_NAME]\_CHECK-PTX-COIL

REFERENCE


MEASUREMENT

MODULE



Reference = saved data from day 1

PHASE



Visual check : no global changes between reference and measurement

# Data comparison

Dicom series [PROTOCOL\_NAME]\_CHECK-PTX-COIL


Predefined ROI

channel by channel comparison

**Correlation (mag+phase)** between reference and measurement in ROI

**RMSE** between reference and measurement in ROI

=> thresholds to be determined





# Data comparison

Dicom series [PROTOCOL\_NAME]\_CHECK

Predefined ROI

channel by channel comparison


**Correlation (mag+phase)** between reference measurement in ROI

**RMSE** between reference and measurement in

=> thresholds to be determined




Follow-up RMSE







Follow-up Corr. Coef. magnitude




Follow-up Corr. Coef. phase



Follow-up linear cor. module




Follow-up linear cor. phase



# Additional info

To reduce variability of B1 map comparison due to phantom positioning relative to isocenter, phantom position is determined using 3 line scans and B1 maps are acquired by shifting the FOV with this data

Phantom profiles in 3 main directions  
And shift estimation using these profiles



# Additional info

In addition to B1 map comparison, a scattering matrix is computed using Dico forward/reflected measurements on every channel. This matrix should inform about transmit coil stability over time. It should also be possible to use the scattering matrix computed by the Siemens calibration (VE12U).

